

**Scharlab S.L.**

Crtra.Polinya a Sentmenat, Km 8,2

08181 Sentmenat

☎ +34937456400

✉ helpdesk@scharlab.com

**CERTIFICATE OF ANALYSIS****AL0162\_21564103/1**

Product Benzyl alcohol, Pharmpur®, Ph Eur, BP, NF

Batch 21564103

**AL0162**

Quality release date 03/03/2021

Expiry date 03/2023

Analysis	Batch value	Specifications
Assay (G.C.)	99,99 %	98,0 - 100,5
Identification	passes test	passes test
Appearance of solution	clear and colourless	clear and colourless
Refractive index n <sub>20</sub> /D	1,5409	1,538 - 1,541
Acidity	passes test	passes test
Peroxide index	0,01	≤ 5
Related substances:		
Impurity A	0,0040 %	≤ 0,15
Impurity B	< 0,0001 %	≤ 0,10
Sum of other peaks	0,0011 %	≤ 0,04
Sum of peaks	0,0134 %	≤ 0,3
Disregard limit	0,0001 %	≤ 0,0001
Residue on evaporation	0,0018 %	≤ 0,05

This product is intended to be suitable for the manufacture of pharmaceuticals, medical devices, cosmetics and food products, but is entirely under the user's responsibility.