


CATALOG NUMBER.....:		002425-0500	
SIZE (NOM VOL mL/oz)....:		500/16	
NOM BRIM CAP (nL).....:		535	
MATERIAL	BOTTLE.....:	LOPE	
	CLOSURE.....:	PP	
NECK FINISH (mm).....:		28-415	
THREAD (T.P.I.).....:		6	
NOM. WT W/CLOSURE (gms).....:		60	
		INCH	MM
A (NECK I.D.).....:	0.84 ± 0.02	21.3 ± 0.5	
H (NOM OVERALL HT)....:	8.3	211	
C (HT. W/O CLOSURE)....:	6.59 ± 0.09	167.4 ± 2.3	
D (DIA).....:	2.86 ± 0.06	72.6 ± 1.5	
CLOSURE COLOR:		CHEMICAL:	
RED		ACETONE	
WHITE		ETHANOL	
BLUE		METHANOL	
YELLOW		ISOPROPANOL	
NATURAL		DISTILLED WATER	

NOTE EACH ASSORTMENT PACK CONSISTS OF 5 BOTTLES, 1 OF EACH OF 5 DIFFERENT WASH BOTTLES WITH RESPECTIVE CHEMICAL NAMES LISTED IN CHART ABOVE.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	NS	APP'D	PVC	PSS 00-2425-0500	REV.	0502
CK'D	PVC		5/6/02			

DESCRIPTION: RIGHT-TO-KNOW
SAFETY WASH BOTTLES ASSORTMENT