


CATALOG NO.	NOMINAL CAPACITY (mL)		MATERIALS		
			VIAL	CLOSURE	GASKET
5000	1.0	1.5	PP	PP	WHITE SILICONE

SIZE CODE.....	-1012		-1020	
NOM. CAP. (mL).....	1.0		1.5	
NOM. BRIM CAP. (mL):	1.5		2.0	
THREAD (T.P.I.).....	12		12	
NOM. WEIGHT W/CLOSURE (gn).....	1.4		1.6	
CLOSURE NOM WEIGHT (gn).....	0.5		0.5	
	in.	mm	in.	mm
A (VIAL NECK I.D.):	0.33 ± 0.01	8.4 ± 0.3	0.33 ± 0.01	8.4 ± 0.3
B (HT. W/CLOSURE):	1.52 ± 0.02	38.6 ± 0.5	1.89 ± 0.02	48.0 ± 0.5
C (HT. W/O CLOSURE):	1.31 ± 0.02	33.3 ± 0.5	1.68 ± 0.02	42.7 ± 0.5
D (VIAL O.D.):	0.42 ± 0.01	10.7 ± 0.3	0.42 ± 0.01	10.7 ± 0.3
F (CLOSURE HT.):	0.50 ± 0.01	12.7 ± 0.3	0.50 ± 0.01	12.7 ± 0.3
G (CLOSURE O.D.):	0.48 ± 0.01	12.2 ± 0.3	0.48 ± 0.01	12.2 ± 0.3
W.T.(NOM. WALL THK):	0.045	1.14	0.045	1.14

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-5000-10XX-1	REV.	0195
CK'D	GRP		1/30/95			

DESCRIPTION:

SYSTEM 100 CRYOGENIC VIALS