


CATALOG NO.	NOMINAL CAPACITY (mL)		MATERIALS	
			VIAL	CLOSURE
005012	1.2	2.0	PP	HDPE


SIZE CODE.....:	-0012		-0020	
SIZE (NOM VOL mL)...:	1.2		2.0	
NOM. BRIM CAP. (mL):	2.0		2.8	
THREAD (T.P.I.).....:	12		12	
NOM. WEIGHT W/CLOSURE (gms).....:	1.5		1.8	
CLOSURE NOM WEIGHT (gms).....:	0.6		0.6	
	IN.	MM	IN.	MM
A (VIAL NECK I.D.)...:	0.39 ±0.01	9.9 ±0.3	0.39 ±0.01	9.9 ±0.3
B (HT. W/CLOSURE)...:	1.51 ±0.02	38.4 ±0.5	1.91 ±0.02	48.5 ±0.5
C (HT. W/O CLOSURE):	1.31 ±0.02	33.3 ±0.5	1.71 ±0.02	43.4 ±0.5
D (VIAL O.D.).....:	0.46 ±0.01	11.7 ±0.3	0.46 ±0.01	11.7 ±0.3
F (CLOSURE HT).....:	0.60 ±0.01	15.2 ±0.3	0.60 ±0.01	15.2 ±0.3
G (CLOSURE O.D.)...:	0.53 ±0.01	13.5 ±0.3	0.53 ±0.01	13.5 ±0.3
W.T.(NOM. WALL THK):	0.040	1.02	0.040	1.02

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CLOSURE DETAILS



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC ²	APP'D	PSS 00-5012-00XX	REV. 0404
CK'D			
DESCRIPTION:		BULK-PACKED STERILE CRYOGENIC VIALS	