


CATALOG NO.....:	5970-0X13		5970-0X16		5970-0X20	
TEST TUBE SIZE (mm):	13		16		20	
MATERIAL.....:	RESMER		RESMER		RESMER	
NOM. WEIGHT (gns):	120		200		190	
	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....:	7.90±0.07	200.7±1.8	9.80±0.07	248.9±1.8	9.90±0.07	251.5±1.8
B (WIDTH).....:	4.06±0.05	103.1±1.3	4.98±0.05	126.5±1.3	4.06±0.05	103.1±1.3
C (HEIGHT).....:	2.30±0.05	58.4±1.3	2.84±0.05	72.1±1.3	3.40±0.05	86.4±1.3
D (SQUARE).....:	0.54±0.01	13.7±0.3	0.69±0.01	17.5±0.3	0.84±0.01	21.3±0.3
HOLE ARRAY.....:	6 x 12		6 x 12		4 x 10	

CATALOG NO.....:	5970-0X25		5970-0X30	
TEST TUBE SIZE (mm):	25		30	
MATERIAL.....:	RESMER		RESMER	
NOM. WEIGHT (gns):	250		220	
	INCH	MM	INCH	MM
A (LENGTH).....:	11.80±0.07	299.7±1.8	11.20±0.07	284.5±1.8
B (WIDTH).....:	4.84±0.05	122.9±1.3	3.38±0.05	85.9±1.3
C (HEIGHT).....:	3.75±0.05	95.3±1.3	4.32±0.05	109.7±1.3
D (SQUARE).....:	1.03±0.01	26.2±0.3	1.25±0.01	31.8±0.3
HOLE ARRAY.....:	4 x 10		3 x 8	

COLOR	CATALOG NUMBER				
WHITE	5970-0013	5970-0016	5970-0020	5970-0025	5970-0030
ORANGE	5970-0113	5970-0116	5970-0120	5970-0125	5970-0130
YELLOW	5970-0213	5970-0216	5970-0220	5970-0225	5970-0230
BLUE	5970-0313	5970-0316	5970-0320	5970-0325	5970-0330
GREEN	5970-0413	5970-0416	5970-0420	5970-0425	5970-0430
RED	5970-0513	5970-0516	5970-0520	5970-0525	5970-0530

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	PVC	PSS 00-5970-0XXX	REV.	0204
CK'D	GRP	2/26/04				
DESCRIPTION: UNWIRE TEST TUBE RACKS						