


CATALOG NUMBER.....:		153732	
NOM. CAP (mL).....:		260	
CAP TYPE.....:		VENT/CLOSE	
CULTURE AREA (cm ²).....:		80	
MATERIAL	FLASK.....:	PS	
	CAP.....:	BLUE HDPE	
		INCH	MM
A (NECK I.D.).....:		0.77	19.7
B (LENGTH W/O CAP).....:		5.62	142.7
C (LENGTH W/ CAP).....:		5.72	145.3
D (BODY WIDTH).....:		3.08	78.2
E (BODY DEPTH).....:		1.44	36.6



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC2	APP'D	PVC	PSS 153732	REV.	0505
CK'D	RDR	6/29/05				

DESCRIPTION: 80cm² NUNCLON FLASK, STERILE