


CATALOG NUMBER.....:		156502	
NOM. CAP. (mL).....:		800	
NOM. TYPE.....:		VENT/CLOSE	
CULTURE AREA (cm <sup>2</sup> ).....:		175	
MATERIAL	FLASK.....:	PS	
	CAP.....:	BLUE HDPE	
		INCH	MM
A (NECK I.D.).....:		0.99	25.1
B (LENGTH W/O CAP).....:		7.64	194.1
C (LENGTH W/ CAP).....:		7.59	192.8
D (BODY WIDTH).....:		4.75	120.7
E (BODY DEPTH).....:		2.01	51.1



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 156502	REV.	0505
CK'D	RDR	6/29/05				
DESCRIPTION: 175cm <sup>2</sup> NUNCLON FLASK, STERILE						