


CARBOY MATERIAL	PC	
CLOSURE MATERIAL	WHITE PP	
SIZE (NOMINAL VOLUME L/GAL)	10/2.5	
NOMINAL BRIM CAPACITY	12.9	
NECK FINISH (mm)	83B	
THREAD TPI	2	
NOMINAL WEIGHT W/CLOSURE (grams)	800	
	INCHES	MM (REF)
A (NECK ID)	2.75±.2	69.9
B (HEIGHT WITH CLOSURE)	15.3±.3	389
C (HEIGHT WITHOUT CLOSURE)	14.9±.2	379
D (DIAMETER)	9.9±.2	251
E (NOMINAL WALL THICKNESS)	.07	1.8



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D PWC	PSS 2251-0020	REV. 0317
CK'D JC	3/2/17		

DESCRIPTION: PELYCHATY

10 LITER ROUND CARBOY