


CATALOG NUMBER.....	SEE BELOW	
NOM. CAP. (µL).....	450	
MATERIAL.....	PP	
	INCH	MM
A (BASE LENGTH).....	5.030	127.76
B (BASE WIDTH).....	3.365	85.47
C (TOP LENGTH).....	4.85	123.3
D (TOP WIDTH).....	3.20	81.2
E (FLANGE LONG SIDE).....	0.30	7.5
F (OVERALL HEIGHT).....	0.57	14.4
G DIA.....	0.40	8.5
H (A-1 LOCATION) (X).....	0.57	14.4
J (A-1 LOCATION) (Y).....	0.44	11.2
K (WELL-TO-WELL SPACING).....	0.35	9.0
L (WELL DEPTH).....	0.43	11.0
M (BOTTOM OF WELL DISTANCE).....	.10	2.5
N (FLANGE SHORT SIDE).....	.09	2.3
O (WELL DIA, BTM OF CYL).....	.30	7.6
P (WELL CONICAL DEPTH).....	.12	3.0

CATALOG NUMBER	COLOR	STERILE	UNITS PER PACK/CASE
249944	NATURAL	NO	20/120
249946	NATURAL	YES	20/120
249945	BLACK	NO	20/120
249949	WHITE	NO	20/120
249943	RED	NO	20/120
249947	YELLOW	NO	20/120
249950	BLUE	NO	20/120

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JDD	APP'D	PVC	PSS 2499XX	REV. 0505
CK'D	RDR 6/29/05		6/29/05		
DESCRIPTION: 96-WELL, PP, "V-BOTTOM" PLATE (NON-TREATED)					