

CATALOG NUMBER.....	SEE BELOW	
NOM. CAP. (µL).....	120	
MATERIAL.....	PS	
	INCH	MM
A (BASE LENGTH).....	5.03	127.8
B (BASE WIDTH).....	3.37	85.5
C (TOP LENGTH).....	4.87	123.7
D (TOP WIDTH).....	3.20	81.3
E (FLANGE LONG SIDE).....	0.25	6.4
F (OVERALL HEIGHT).....	0.57	14.4
G DIA.....	0.15	3.7
H (A-1 LOCATION) (X).....	0.48	12.1
J (A-1 LOCATION) (Y).....	0.35	9.0
K (WELL-TO-WELL SPACING)...	0.18	4.5
L (WELL DEPTH).....	0.46	11.7
M (BOTTOM OF WELL DISTANCE):	.07	1.7
N (FLANGE SHORT SIDE).....	.10	2.5
*P (LID LENGTH).....	5.02	127.4
*Q (LID WIDTH).....	3.36	85.3
*R (LID HEIGHT).....	.36	9.1
S (STACKED HEIGHT - NO LID):	1.09	27.7
*T (STACKED HEIGHT W/LIDS):	1.33	33.8
*U (ASSEMBLY HEIGHT W/LID):	.665	16.9




CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
164564	CELL CULTURE	BLACK	YES	YES	10/30
164610	CELL CULTURE	WHITE	YES	YES	10/30
164688	CELL CULTURE	CLEAR	YES	YES	10/30
164571	CELL CULTURE	BLACK	YES	NO	25/100
165195	CELL CULTURE	WHITE	YES	NO	25/100
164555	CELL CULTURE	CLEAR	YES	NO	25/100
264556	NON-TREATED	BLACK	NO	NO	10/30
264572	NON-TREATED	WHITE	NO	NO	10/30
265196	NON-TREATED	CLEAR	NO	NO	10/30
242757	NON-TREATED	CLEAR	YES	YES	10/30

CATALOG NO.	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
242765	NON-TREATED	CLEAR	YES	YES	1/30
262260	NON-TREATED	BLACK	NO	NO	25/100
262360	NON-TREATED	WHITE	NO	NO	25/100
262160	NON-TREATED	CLEAR	NO	NO	25/100
265202	NON-TREATED	CLEAR	YES	YES	25/100
265203	NON-TREATED	CLEAR	YES	NO	25/100
460435	POLYSORP®	BLACK	NO	NO	10/30
460518	MAXISORP™	BLACK	NO	NO	10/30
460372	MAXISORP™	WHITE	NO	NO	10/30
464718	MAXISORP™	CLEAR	NO	NO	10/30

* DIM.S APPLY ONLY WHERE "WITH LID" = YES, USING NUNC LID #8337.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JDD	APP'D	GRP	PSS 1646XX	REV. 0706
CK'D	AEG		7/24/06		

DESCRIPTION:
384 WELL, PS, FLAT BOTTOM PLATE