


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (μL)	50-1700μL	
MATERIAL	PP	
	MM	INCH
A, BASE LENGTH	127,8	5,03
B, BASE WIDTH	85,5	3,37
C, TOP LENGTH	123,6	4,87
D, TOP WIDTH	81,6	3,21
E, FLANGE LONG SIDE	2,5	0,1
F, OVERALL HEIGHT	43,6	1,72
G1, DIA TOP	Ø8,5	Ø0,33
G2, DIA BOTTOM	Ø7,2	Ø0,28
H, A - I LOCATION (X)	14,3	0,56
J, A - I LOCATION (Y)	11,3	0,44
K, WELL TO WELL SPACING	9	0,35
L, WELL DEPTH	41,5	1,63
M, BOTTOM OF WELL DISTANCE	2,5	0,1
N, FLANGE SHORT SIDE	2,5	0,1

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/ PK/CASE
278743	NON-TREATED	NATURAL	+	-	1/60
278752	NON-TREATED	NATURAL	-	-	5/60

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
DK-4000 Roskilde • Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT: U96 DEEPWELL PLATE 2.0ML., PP		DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORMAT:
ITEM NUMBER:	MATERIAL:	DRAWINGS NO.:	VERSION:	PAGE TYPE:	SCALE:	METRIC Metric Preference

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC