


CATALOG NUMBER	SEE BELOW	
WORKING RANGE (µL)	10 – 250 µL	
MATERIAL	PP	
	MM	INCH
A , BASE LENGTH	125,9	5
B , BASE WIDTH	84,7	3,33
C , TOP LENGTH	122,8	4,83
D , TOP WIDTH	81,1	3,19
E , FLANGE LONG SIDE	2,7	0,11
F , OVERALL HEIGHT	14,4	0,57
G1 , DIA TOP	Ø6,7	Ø0,26
G2 , DIA BOTTOM	Ø6,4	Ø0,25
H , A - I LOCATION (X)	14,2	0,56
J , A - I LOCATION (Y)	11,3	0,45
K , WELL TO WELL SPACING	8,9	0,35
L , WELL DEPTH	9,8	0,39
M , BOTTOM OF WELL DISTANCE	4,6	0,18
N , FLANGE SHORT SIDE	2,7	0,11

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	UNITS/PK/CASE
442587	NON-TREATED	NATURAL	-	-	5/60

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



ThermoFisher
SCIENTIFIC

Roskilde Site · Kamstrupvej 90 · P.O.Box 280
DK-4000 Roskilde · Denmark
Tel +45 4631 2000 · Fax +45 4631 2042
E-mail: info.nunc@thermofisher.com

www.nuncbrand.com

SUBJECT: V96 MICROWELL PLATE, PP	DRAWN BY/DATE:	CHANGE ORDER NO.:	VALID FROM:	APPROVED BY/DATE:	FORM/
ITEM NUMBER:	MATERIAL:	DRAWING NO.:	VERSION:	PAGE TYPE:	SCALE:

CONFIDENTIAL
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC