

NOT TO SCALE

127.76 ±0.25

99.00 ±0.25

14.38 ±0.25

8.40 ±0.10

(x96)

9.00

Y

Copyright ABgene®

This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.

21.00 ±0.10

0.50



SECTION Y-Y

8.40 ±0.10 (x96)

Y

46.00

PLATES STACKED



SECTION X-X

Pub No. MAN0014419

All dimensions in mm

TITLE		
1.2ml Square Well Storage Plate		
DRG. NO.	ISSUE	DATE
AB-1127	D	20/9/02