

NOT TO SCALE

Pub. no. MAN0014431

127.76 ±0.25

99.00 ±0.25

14.38 ±0.25

9.00

Y

Ø7.00 ±0.10
(x96)

39.35 ±0.10

0.50

4.00 (X8)

25.00

25.00

85.48 ±0.25


63.00 ±0.25

11.24 ±0.25


9.00

X

X



SECTION Y-Y




34.00

34.00

Y


42.50 ±0.25

3.00 ±0.25



SECTION X-X

PLATES STACKED



TITLE

1.2ml Storage Plate

DRG. NO.

AB-0564

ISSUE

E

DATE

7/5/02

All dimensions in mm

Copyright ABgene®

This drawing and/or any technical information relating to it is issued in strict confidence and may not be reproduced or used for manufacture without written permission. Upon receipt of which it is the users responsibility to seek confirmation of any dimensions and production tolerances.